

Lapiplasty®

3D Bunion Correction™


Bunions are more than a bump.

A common misconception is that a bunion is simply an overgrowth of bone to be "shaved off".

However, a bunion is a complex, 3D deformity caused by an unstable joint in the middle of the foot which allows the metatarsal bone to drift out of alignment and to **1 lean, 2 elevate & 3 rotate.**

Bunion "bump" appears here

Root cause of the bunion - an unstable joint - is actually here



Traditional Osteotomy Surgery vs. Lapiplasty® 3D Bunion Correction™

- X Unnaturally cuts & shifts bone; only a 2D correction
- X Addresses cosmetic "bump" only; not the root cause
- X Patients may be off their feet for up to 6-8 weeks
- X High recurrence; up to 70% of bunions can return¹

- ✓ Returns entire bone to normal alignment; a 3D correction
- ✓ Fixes the root cause of the bunion; an unstable joint
- ✓ Get back on your feet quickly; in many cases within 2 weeks²
- ✓ Low recurrence; 97% maintain 3D correction @ 13 months³


How it works

Traditional 2D Osteotomy

Shaves the bump, then cuts & shifts top part of bone over. This addresses the cosmetic bump, but does not address the root cause of the bunion; an unstable joint. The bone is now crooked with the bunion likely to return.¹

1


Shaves off bump



Unstable Foundation
(root of the problem is not corrected)

2

Cuts & shifts top of bone




Unstable Foundation
(root of the problem is not corrected)

Lapiplasty® 3D Bunion Correction™

Rather than cutting a straight metatarsal bone, the entire bone is returned to natural, anatomic alignment in all 3 dimensions. Patented titanium plates permanently secure the unstable joint and allow early weight-bearing.

1


Returns entire bone to normal 3D alignment



3D Bunion Correction™
(normal position in all 3 dimensions)

2


Permanently secures unstable joint



Joint Secured
(unstable joint permanently secured)



Post-Op - Interior
X Crooked Bone, Unstable Joint




Post-Op - Exterior
X Only Pretty On Outside

up to **70%**
osteotomy patients
have their bunions
return over time¹

Clinical Data
X High Recurrence¹



Post-Op - Interior
✓ Straight Bone, Joint Secured



Post-Op - Exterior
✓ Pretty Inside & Out

97%
Lapiplasty® patients
retain their 3D
correction³
in a 13 month study

Clinical Data
✓ Low Recurrence³

As with any medical treatment, individual results may vary and these experiences are unique and specific to these patients. And, not every patient is a Lapiplasty® candidate.

Only a surgeon can tell if the Lapiplasty® Procedure is right for you. There are potential risks and recovery takes time.

Potential risks include but are not limited to infection, discomfort from the presence of the implant, loosening of the implant, and loss of correction with nonunion or malunion.

Lapiplasty® is a registered mark of Treace Medical Concepts, Inc. Fix it Right the First Time™ and 3D Bunion Correction™ are trademarks of Treace Medical Concepts, Inc. Copyright 2019 Treace Medical Concepts, Inc. All rights reserved. M052A

1 Jeuken RM, et al. Foot Ankle Int. 2016; 37:687-95. 2 Smith B, et al. 2017 AOFAS Annual Meeting, Seattle, WA. 3 Santrock R, et al. 2018 AOFAS Annual Meeting, Boston, MA.